

Rabbit Anti-Mouse IgG (H&L) - AF790

Catalog No: #L35639

Package Size: #L35639-1 100 ul #L35639-2 500 ul #L35639-3 1 ml

Description

Product Name	Rabbit Anti-Mouse IgG (H&L) - AF790
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	M
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Mouse IgG
Conjugates	AF Dye
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.