

Goat Anti-Human IgA - HRP

Catalog No: #L35662

Package Size: #L35662-1 100ul #L35662-2 500ul #L35662-3 1ml

Description

Product Name	Goat Anti-Human IgA - HRP
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	H
Immunogen Description	Human IgA
Conjugates	HRP
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Proclin
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.