

Goat Anti-Mouse IgG2a - AF488

Catalog No: #L35675

Package Size: #L35675-1 100ul

Description

Product Name	Goat Anti-Mouse IgG2a - AF488
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	ELISA;IHC;IF;FC
Species Reactivity	Mouse
Immunogen Description	Mouse IgG2a
Conjugates	AF488 Ex/Em = 490/515 nm
Concentration	1 mg/ml
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

IHC 1:100-1:1000;IF 1:100-1:1000;FC 1:1000-1:4000;ELISA Use at an assay dependent concentration

Note: This product is for in vitro research use only and is not intended for use in humans or animals.