

UMPS (Phospho-Tyr37) Antibody

Catalog No: #SAB319P

Package Size: #SAB319P-1 50ul #SAB319P-2 100ul

Description

Product Name	UMPS (Phospho-Tyr37) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	Custom antibody
Species Reactivity	Hu
Immunogen Type	Peptide-KLH
Conjugates	Unconjugated
Target Name	UMPS
Other Names	Uridine 5'-monophosphate synthase; UMP synthase
Accession No.	uniprot:P11172
Calculated MW	52kDa
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Western blotting: 1:500~1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.