

PAQR4 Antibody

Catalog No: #SAB607

Package Size: #SAB607 100ul

Description

Product Name	PAQR4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified from rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns.
Applications	WB
Species Reactivity	Human,Mouse
Specificity	PAQR4 Antibody detects endogenous levels of PAQR4.
Immunogen Description	A synthesized peptide derived from human PAQR4.
Other Names	Progestin and adipoQ receptor family member 4,Progestin and adipoQ receptor family member IV,
SDS-PAGE MW	29kDa
Concentration	1 mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl,0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000

Background

Molecular function for PAQR4 involves signaling receptor activity. "Phylogenetic-based propagation of functional annotations within the Gene Ontology consortium."

Gaudet P., Livstone M.S., Lewis S.E., Thomas P.D.

Brief Bioinform 12:449-462(2011)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.