

Recombinant HA-Tag Monoclonal Antibody(12H2),
HRP Conjugated

Catalog No: #T403

Package Size: #T403-1 50ul #T403-2 100ul



Description

Product Name	Recombinant HA-Tag Monoclonal Antibody(12H2), HRP Conjugated
Clonality	Monoclonal
Purification	affinity purified by Protein A
Applications	WB,ELISA
Species Reactivity	Species independent
Conjugates	HRP
Formulation	1Γ PBS (pH7.4) αO 0.05% proclin300αO 50% Glycerol
Storage	Store at -20°C/1 year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:30000-1:100000, ELISA 1:5000-1:10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.