

Recombinant His-Tag Monoclonal Antibody(14W7),  
HRP Conjugated

Catalog No: #T404

Package Size: #T404-1 50ul #T404-2 100ul



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

|                    |                                                               |
|--------------------|---------------------------------------------------------------|
| Product Name       | Recombinant His-Tag Monoclonal Antibody(14W7), HRP Conjugated |
| Clonality          | Monoclonal                                                    |
| Purification       | affinity purified by Protein A                                |
| Applications       | WB,ELISA                                                      |
| Species Reactivity | Species independent                                           |
| Conjugates         | HRP                                                           |
| Formulation        | 1Γ PBS (pH7.4) αO 0.05% proclin300αO 50% Glycerol             |
| Storage            | Store at -20°C/1 year. Avoid repeated freeze/thaw cycles.     |

Application Details

WB 1:30000-1:100000, ELISA 1:5000-1:10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.