

Strep II tag Mouse Monoclonal Antibody

Catalog No: #T603

Package Size: #T603 100ul

Description

Product Name	Strep II tag Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Affinity-chromatography
Applications	WB
Specificity	Strep II tag mouse monoclonal antibody detects endogenous levels of internal, C-terminal, and N-terminal Strep II tagged proteins.
Immunogen Type	Peptide
Immunogen Description	A synthetic peptide NWSHPQFEK.
Conjugates	Unconjugated
Target Name	Strep II tag
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20 °C. Stable for 12 months from date of receipt

Application Details

WB: 1:3000~1:10000

IF/ICC: 1:200

Images



1ug Strep II fusion protein+ Primary antibody dilution at
1 μ g • 1:5,000 2 μ g • 1:10,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.