

AmCyan Monoclonal Antibody

Catalog No: #T614

Package Size: #T614-1 50ul #T614-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	AmCyan Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Specificity	The antibody detects AmCyan & AmCyan Tag protein.
Conjugates	Unconjugated
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:5000

Images



Western blot analysis of recombinant AmCyan protein, diluted at 1:5000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.